

Containers Container



Info

Suitability of our sterile containers for steam sterilization procedures¹⁾ as per DIN 58953 part 9 (issue Mai 1987)

Sterilisation container	Packing		Steam sterilization procedure acc. to DIN 58946, part I				
	Description	Syboic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injektion procedure
with filters (StB-F)	Textiles in sterilation container with filter in lid and bottom acc. to DIN 58952 part I		+	+	+	+	+
			○	○	+	+	+
	Textiles in sterilation container with filter in lid, bottom not perforated acc. to DIN 58952 part I		-	-	+	+	+
			-	-	+	+	○
	Textiles in sterilation bottom, lid not perforated acc. to DIN 58952 part I		○	+	+	+	+
			-	○	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III in sterilation container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc. to DIN 58952 part III, in sterilation container with filter in lid and bottom not perforated acc. to DIN 58952 part I		-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc. to DIN 58952 part III, in sterilation container with filter in bottom acc. to DIN 58952 part I		○	+	+	+	+
	with Valves (StB-V)	Textiles in sterilation container with valves in lid and bottom acc. to DIN 58952 part I		-	○	+	+
			-	-	+	+	○
Textiles in sterilation container with valves in lid no valves on the bottom			-	-	+	+	-
			-	-	+	+	-
Instruments in wire basket in sterilization container with valve in lid and bottom			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
Instruments in wire basket in sterilization container with valve in lid and not perforated			-	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols: + applicable - not applicable ○ limited application

¹⁾ In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.

²⁾ Suitable are only sterilization containers guaranteeing evaporation or discharge of condensate.

³⁾ Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.

Grosse Container und Zubehör
Large Containers and Accessories



Grosse Container und Zubehör

Large Containers and Accessories

Basismodell Basic Model


Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated


Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



580 x 280 x 100 mm








-  ER110.100
-  ER120.100
-  ER130.100
-  ER140.100
-  ER150.100
-  ER160.100

-  ER111.100
-  ER121.100
-  ER131.100
-  ER141.100
-  ER151.100
-  ER161.100

580 x 280 x 135 mm



-  ER110.135
-  ER120.135
-  ER130.135
-  ER140.135
-  ER150.135
-  ER160.135

-  ER111.135
-  ER121.135
-  ER131.135
-  ER141.135
-  ER151.135
-  ER161.135

580 x 280 x 150 mm



-  ER110.150
-  ER120.150
-  ER130.150
-  ER140.150
-  ER150.150
-  ER160.150

-  ER111.150
-  ER121.150
-  ER131.150
-  ER141.150
-  ER151.150
-  ER161.150

580 x 280 x 200 mm







-  ER110.200
-  ER120.200
-  ER130.200
-  ER140.200
-  ER150.200
-  ER160.200

-  ER111.200
-  ER121.200
-  ER131.200
-  ER141.200
-  ER151.200
-  ER161.200

580 x 280 x 260 mm



-  ER110.260
-  ER120.260
-  ER130.260
-  ER140.260
-  ER150.260
-  ER160.260

-  ER111.200
-  ER121.200
-  ER131.200
-  ER141.200
-  ER151.200
-  ER161.200



Grosse Container und Zubehör

Large Containers and Accessories

Sicherheitsdeckel Safe Model



Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



580 x 280 x 100 mm



-  ER112.100
-  ER122.100
-  ER132.100
-  ER142.100
-  ER152.100
-  ER162.100

-  ER113.100
-  ER123.100
-  ER133.100
-  ER143.100
-  ER153.100
-  ER163.100

580 x 280 x 135 mm



-  ER112.135
-  ER122.135
-  ER132.135
-  ER142.135
-  ER152.135
-  ER162.135

-  ER113.135
-  ER123.135
-  ER133.135
-  ER143.135
-  ER153.135
-  ER163.135

580 x 280 x 150 mm



-  ER112.150
-  ER122.150
-  ER132.150
-  ER142.150
-  ER152.150
-  ER162.150

-  ER113.150
-  ER123.150
-  ER133.150
-  ER143.150
-  ER153.150
-  ER163.150

580 x 280 x 200 mm








-  ER112.200
-  ER122.200
-  ER132.200
-  ER142.200
-  ER152.200
-  ER162.200

-  ER113.200
-  ER123.200
-  ER133.200
-  ER143.200
-  ER153.200
-  ER163.200

580 x 280 x 260 mm



-  ER112.260
-  ER122.260
-  ER132.260
-  ER142.260
-  ER152.260
-  ER162.260

-  ER113.260
-  ER123.260
-  ER133.260
-  ER143.260
-  ER153.260
-  ER163.260

Grosse Container und Zubehör

Large Containers and Accessories

Zubehör

Accessories

Siebkörbe
Wire Baskets




485 x 255 x 30 mm **ER190.030**
485 x 255 x 50 mm **ER190.050**
485 x 255 x 70 mm **ER190.070**
485 x 255 x 100 mm **ER190.100**

Siebkörbe
Wire Baskets



540 x 255 x 30 mm **ER191.030**
540 x 255 x 50 mm **ER191.050**
540 x 255 x 70 mm **ER191.070**
540 x 255 x 100 mm **ER191.100**

Siebboden
Wire Base



545 x 255 x 30 mm **ER193.000**

Silicon
Silicon Mat




445 x 235 cm **ER852.000**

Spannrahmen
Drape Retainer



545 x 255 mm **ER195.000**

Tuch
Drape



130 x 60 cm **ER095.070**
140 x 100 cm **ER095.080**
185 x 45 cm **ER095.090**
185 x 60 cm **ER095.100**

3/4 Container und Zubehör
3/4 Containers and Accessories



3/4 Container und Zubehör

3/4 Containers and Accessories

Basismodell Basic Model

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



465 x 280 x 100 mm











-  ER210.100
-  ER220.100
-  ER230.100
-  ER240.100
-  ER250.100
-  ER260.100

-  ER211.100
-  ER221.100
-  ER231.100
-  ER241.100
-  ER251.100
-  ER261.100

465 x 280 x 135 mm








-  ER210.135
-  ER220.135
-  ER230.135
-  ER240.135
-  ER250.135
-  ER260.135

-  ER211.135
-  ER221.135
-  ER231.135
-  ER241.135
-  ER251.135
-  ER261.135

465 x 280 x 150 mm



-  ER210.150
-  ER220.150
-  ER230.150
-  ER240.150
-  ER250.150
-  ER260.150

-  ER211.150
-  ER221.150
-  ER231.150
-  ER241.150
-  ER251.150
-  ER261.150



3/4 Container und Zubehör

3/4 Containers and Accessories

Sicherheitsdeckel Safe Model

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



465 x 280 x 100 mm






-  ER212.100
-  ER222.100
-  ER232.100
-  ER242.100
-  ER252.100
-  ER262.100

-  ER213.100
-  ER223.100
-  ER233.100
-  ER243.100
-  ER253.100
-  ER263.100

465 x 280 x 135 mm








-  ER212.135
-  ER222.135
-  ER232.135
-  ER242.135
-  ER252.135
-  ER262.135

-  ER213.135
-  ER223.135
-  ER233.135
-  ER243.135
-  ER253.135
-  ER263.135

465 x 280 x 150 mm



-  ER212.150
-  ER222.150
-  ER232.150
-  ER242.150
-  ER252.150
-  ER262.150

-  ER213.150
-  ER223.150
-  ER233.150
-  ER243.150
-  ER253.150
-  ER263.150

3/4 Container und Zubehör

3/4 Containers and Accessories

Zubehör


Accessories

Siebkörbe
Wire Baskets



405 x 255 x 30 mm **ER290.030**
405 x 255 x 50 mm **ER290.050**
405 x 255 x 70 mm **ER290.070**
405 x 255 x 100 mm **ER290.100**

Siebboden
Wire Base




405 x 255 x 30 mm **ER293.030**

Spannrahmen
Drape Retainer



405 x 255 mm **ER295.000**

Drape
Tuch



140 x 100 cm **ER095.080**

Halb-Container und Zubehör
Half-Size-Containers and Accessories



Halb-Container und Zubehör

Half-Size-Containers and Accessories

Basismodell Basic Model

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



285 x 280 x 100 mm



-  ER310.100
-  ER320.100
-  ER330.100
-  ER340.100
-  ER350.100
-  ER360.100

-  ER311.100
-  ER321.100
-  ER331.100
-  ER341.100
-  ER351.100
-  ER361.100

285 x 280 x 135 mm



-  ER310.135
-  ER320.135
-  ER330.135
-  ER340.135
-  ER350.135
-  ER360.135

-  ER311.135
-  ER321.135
-  ER331.135
-  ER341.135
-  ER351.135
-  ER361.135

285 x 280 x 150 mm



-  ER310.150
-  ER320.150
-  ER330.150
-  ER340.150
-  ER350.150
-  ER360.150

-  ER311.150
-  ER321.150
-  ER331.150
-  ER341.150
-  ER351.150
-  ER361.150

285 x 280 x 200 mm



-  ER310.200
-  ER320.200
-  ER330.200
-  ER340.200
-  ER350.200
-  ER360.200

-  ER311.200
-  ER321.200
-  ER331.200
-  ER341.200
-  ER351.200
-  ER361.200

285 x 280 x 260 mm



-  ER310.260
-  ER320.260
-  ER330.260
-  ER340.260
-  ER350.260
-  ER360.260

-  ER311.260
-  ER321.260
-  ER331.260
-  ER341.260
-  ER351.260
-  ER361.260



Halb-Container und Zubehör

Half-Size-Containers and Accessories

Sicherheitsdeckel Safe Model

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



285 x 280 x 100 mm









-  ER312.100
-  ER322.100
-  ER332.100
-  ER342.100
-  ER352.100
-  ER362.100

-  ER313.100
-  ER323.100
-  ER333.100
-  ER343.100
-  ER353.100
-  ER363.100

285 x 280 x 135 mm



-  ER312.135
-  ER322.135
-  ER332.135
-  ER342.135
-  ER352.135
-  ER362.135

-  ER313.135
-  ER323.135
-  ER333.135
-  ER343.135
-  ER353.135
-  ER363.135

285 x 280 x 150 mm



-  ER312.150
-  ER322.150
-  ER332.150
-  ER342.150
-  ER352.150
-  ER362.150

-  ER313.150
-  ER323.150
-  ER333.150
-  ER343.150
-  ER353.150
-  ER363.150

285 x 280 x 200 mm



-  ER312.200
-  ER322.200
-  ER332.200
-  ER342.200
-  ER352.200
-  ER362.200

-  ER313.200
-  ER323.200
-  ER333.200
-  ER343.200
-  ER353.200
-  ER363.200

285 x 280 x 260 mm



-  ER312.260
-  ER322.260
-  ER332.260
-  ER342.260
-  ER352.260
-  ER362.260

-  ER313.260
-  ER323.260
-  ER333.260
-  ER343.260
-  ER353.260
-  ER363.260

Halb-Container und Zubehör

Half-Size-Containers and Accessories

Zubehör

Accessories

Siebkörbe
Wire Baskets



255 x 245 x 30 mm **ER390.030**
255 x 245 x 50 mm **ER390.050**
255 x 245 x 70 mm **ER390.070**
255 x 245 x 100 mm **ER390.100**

Siebboden
Wire Base




258 x 258 x 30 mm **ER393.000**

Spannrahmen
Drape Retainer



255 x 255 mm **ER395.000**

Tuch
Drape



90 x 60 cm **ER095.040**
90 x 80 cm **ER095.050**
130 x 45 cm **ER095.060**
130 x 60 cm **ER095.070**

Breit-Container und Zubehör

Wide Body Containers and Accessories

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



600 x 400 x 120 mm



-  ER110.62
-  ER120.62
-  ER130.62
-  ER140.62
-  ER150.62
-  ER160.62

-  ER111.62
-  ER121.62
-  ER131.62
-  ER141.62
-  ER151.62
-  ER161.62

600 x 400 x 180 mm



-  ER110.68
-  ER120.68
-  ER130.68
-  ER140.68
-  ER150.68
-  ER160.68

-  ER111.68
-  ER121.68
-  ER131.68
-  ER141.68
-  ER151.68
-  ER161.68

Zubehör

Accessories

Siebkorb
Wire Basket




470 X 337 x 51 mm **ER190.55**

Siebkorb
Wire Basket



533 X 337 x 101 mm **ER191.60**


Weitere Container Other Containers



Flach-Container und Zubehör Flat Containers and Accessories

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



285 x 280 x 55 mm



285 x 280 x 85 mm



-  ER310.050
-  ER320.050
-  ER330.050
-  ER340.050
-  ER350.050

-  ER310.080
-  ER320.080
-  ER330.080
-  ER340.080
-  ER350.080

-  ER311.080
-  ER321.080
-  ER331.080
-  ER341.080
-  ER351.080

Zubehör Accessories

Siebkorb
Wire Baskets




255 x 245 x 30 mm **ER390.030**
255 x 245 x 50 mm **ER390.050**

Siebboden
Wire Base



245 x 255 x 30 mm **ER393.000**

Silicon
Silicon Mat



245 x 245 mm **ER851.000**

Spannrahmen
Drape Retainer



255 x 255 mm **ER395.000**

Tuch
Drape




90 x 60 cm **ER095.040**
90 x 80 cm **ER095.050**

Mini Container und Zubehör

Mini Containers and Accessories

Boden ungelocht / Deckel gelocht
Bottom non-perforated / Lid perforated

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



300 x 140 x 40 mm



-  ER410.040
-  ER420.040
-  ER430.040
-  ER440.040
-  ER450.040

-  ER411.040
-  ER421.040
-  ER431.040
-  ER441.040
-  ER451.040

300 x 140 x 70 mm




-  ER410.070
-  ER420.070
-  ER430.070
-  ER440.070
-  ER450.070

-  ER411.070
-  ER421.070
-  ER431.070
-  ER441.070
-  ER451.070

300 x 140 x 100 mm



-  ER410.100
-  ER420.100
-  ER430.100
-  ER440.100
-  ER450.100

-  ER411.100
-  ER421.100
-  ER431.100
-  ER441.100
-  ER451.100

Zubehör Accessories

Silicon
Silicon Mat



270 x 125 mm **ER850.000**

Spannrahmen
Drape Retainer



280 x 125 mm **ER495.000**

Tuch
Drape



60 x 40 cm **ER095.010**

Spezial Container und Zubehör Special Containers and Accessories

Boden gelocht / Deckel gelocht
Bottom perforated / Lid perforated



310 x 190 x 40 mm



ER510.040

310 x 190 x 65 mm



ER510.060

310 x 190 x 130 mm



ER510.130


Zubehör Accessories

Spannrahmen
Drape Retainer



310 x 190 mm **ER595.000**

Tuch
Drape



60 x 40 cm **ER095.010**
65 x 65 cm **ER095.020**
85 x 85 cm **ER095.030**

Spezial Container Special Containers

Boden gelocht / Deckel gelocht (ohne Filter)
Bottom perforated / Lid perforated (without Filter)



300 x 138 x 65 mm



- ER610.060
- ER620.060
- ER630.060
- ER640.060
- ER650.060

500 x 155 x 75 mm



- ER610.070
- ER620.070
- ER630.070
- ER640.070
- ER650.070

Boden gelocht / Deckel gelocht (ohne Filter)
Bottom perforated / Lid perforated (without Filter)



250 x 60 x 50 mm



- ER600.010


450 x 70 x 70 mm



- ER600.020

Spezial Container Special Containers

Boden ungelocht / Deckel ungelocht
Bottom non-perforated / Lid non-perforated



285 x 280 x 100 mm



285 x 280 x 150 mm



580 x 280 x 150 mm



- ER315.100
- ER325.100
- ER335.100
- ER345.100
- ER355.100
- ER365.100

- ER315.150
- ER325.150
- ER335.150
- ER345.150
- ER355.150
- ER365.150

- ER115.150
- ER125.150
- ER135.150
- ER145.150
- ER155.150
- ER165.150

Zubehör Accessories



Zubehör Accessories

Einweg/ Papierfilter
Disposable paper filters



ER800.000

100 Stück, für Große, 3/4- und Halb Container
100 Pcs/Package for Large, 3/4 and Half- Size- Containers

Textile Filters
Textile Filters



ER800.010

Einzel, für Große, 3/4- und Halb- Container
Single, for Large, 3/4 and Half- Size- Containers

Indikatoretiketten
Paper Labels with Indicator



ER810.000

100 Stück, für Große, 3/4- und Halb- Container
100 Pcs., for Large, 3/4 and Half- Size- Containers



ER802.000

100 Stück, für Mini Container
100 Pcs., for Mini Containers

ER804.000

100 Stück, für Dental Container
100 Pcs., for Dental Containers



ER802.010

Einzel, für Mini Container
Single, for Mini Containers

ER804.010

Einzel, für Dental Container
Single, for Dental Containers

Papierfilterspender
Filter Dispenser



ER801.000

Filterhalterung
Filter Fixation



ER809.000

Für Papierfilter
For Paper Filters

ER809.010

Für Textilfilter
For Textile Filters


Verschlussplombe
Plastic Security Seal



ER819.000

100 Stück
100 Pcs.

Kennzeichnungsschilder
Identification Labels



Für max. 16 Zeichen
For max. 16 letters



Spender für Kennzeichnungsschilder
Dispenser for Identification Labels



ER829.010

Zubehör Accessories



ER830.000

Haltestift, für Siebkörbe
Holding Pin, for Wire Baskets



ER830.010

Haltestift, für Trennbleche
Holding Pin, for Partition Sheets




Trennbleche, rostfreier Stahl
Partition Sheets, Stainless Steel

ER831.050 50 x 20 mm

ER831.130 130 x 20 mm

ER831.220 225 x 20 mm

ER831.460 460 x 20 mm



ER840.010

Gleitschieber, für 60-684.02
Divider, for 60-684.02



ER841.000 232 mm

ER841.010 471 mm

Schiene für Gleitschieber, mit Sicherungsschrauben
Limiting Bar, with Securing Screws



ER842.000 232 mm

ER842.010 471 mm

Anschlagrahmen, mit Sicherungsschrauben
Limiting Bar, with Securing Screws



ER845.000

Spiralfeder Haltevorrichtung
Spiral Holding



ER847.000

Doppel- Spiralfeder Haltevorrichtung
Double Spiral Holding



ER849.000

Fixierklammer
Fixation Clamp



Wagen Trolleys



ER890.000

720 x 620 x 800 mm



ER891.000

1030 x 620 x 1170 mm

Siebkörbe mit Deckel Wire Baskets with Lids



Deckel, mit Verschluss für Siebkorb
Lid, with catch for Wire Baskets

ES390.000	243 x 253 mm
ES290.000	405 x 253 mm
ES190.000	485 x 253 mm
ES191.000	540 x 253 mm

Siebkörbe
Wire Baskets

ES390.030	243 x 253 x 30 mm
ES390.050	243 x 253 x 50 mm
ES390.070	243 x 253 x 70 mm
ES390.100	243 x 253 x 100 mm

ES190.030	485 x 253 x 30 mm
ES190.050	485 x 253 x 50 mm
ES190.070	485 x 253 x 70 mm
ES190.100	485 x 253 x 100 mm

ES290.030	405 x 253 x 30 mm
ES290.050	405 x 253 x 50 mm
ES290.070	405 x 253 x 70 mm
ES290.100	405 x 253 x 100 mm

ES191.030	540 x 253 x 30 mm
ES191.050	540 x 253 x 50 mm
ES191.070	540 x 253 x 70 mm
ES191.100	540 x 253 x 100 mm